

QuadFuge™ Plate Spinner is the centrifuge designed specifically for quick and easy centrifugation of samples in PCR plates or microplates. Two standard microplates are loaded through a slot in the top of the centrifuge unit. Sealed plates are loaded vertically.




Cat. No. : QF-2228

Features :

- The speed is adjustable between 2,200 and 2,800 rpm.
- Short spin function is provided for quick spin down.
- Operation can be timed from 15 seconds to 1 minute and 59 seconds
- Accelerates and stops quickly
- Running safely - the rotors stop once the lid release button is pushed
- Brushless, free maintenance.

Specification :

Speed	2,200 ~ 2,800Rpm
Maximum RCF	550 x g
Timer	15 sec. ~ 1 min. 59 sec.
Maximum capacity	2 standard PCR or 96-well plates
Acceleration	≤10sec. (0 to 2,800 rpm)
Deceleration	4 sec. (2,800 to 0 rpm)
Dimensions	26 x 24 x 22.6 cm (W x D x H)
Weight	3.3 kg
Electrical	AC 110 - 220V / 50 - 60 Hz



Slot → Shifting Speed → Spin