



Features :

Liquid aspiration: the user can choose HIGH or LOW work speed and additionally adjust the suction by the pressure applied to the trigger button.

Dispensing: carried out by gravity or supported by pump which empties the pipette with blow out.

Overfilling protection: supplied with PTFE filter blocking any liquid from entering the unit. To protect the samples from cross contamination, filters, silicone resin and nosepiece can be easily exchanged and autoclaved.

Powerful: environmentally friendly 3 Ni-MH batteries enables around 2,000 times of continuous work. The red light indicates when AgileAid should be recharged.

Each unit comes with charging stand, power supply, wall mounting and PTFE filter

Cat. Number	Description
CPA-100	Comfyaid, universal power supply, charging stand, 2 PTFE autoclavable filter 0,2 µm, wall mounting,
CPA-100-1	Silicone resin
CPA-100-2	PTFE autoclavable filter 0.2 µm, 5 pieces/pack
CPA-100-3	PTFE autoclavable filter 0.45 µm, 5 pieces/pack
CPA-100-4	Wall mounting
CPA-100-5	Pipette holder