

## DNA Marker

### 1 kb Ladder, 13 Fragments

Product contents:

P/N: E3011

1x500ul

Sufficient for approximately 100 loads

---

#### Introduction:

Avans DNA Marker is 13 DNA fragment ladder covering wide range of 250-10,000 bp sizes. Conveniently delivered in ready-to-use format and ideal for agarose gel electrophoresis, Northern, Southern, and related DNA analysis applications.

---

#### Protocol:

Swirl tube well.

Add ~5µl per load onto agarose gel prior to electrophoresis.

---

#### Shipping and Storage Conditions :

Shipped in RT

Store in 4 °C for short term

Store in -20 °C for long term

