

DNA Marker

1 kb Plus Ladder, 19 Fragments

Product contents:

P/N: E3012

1x500ul

Sufficient for approximately 100 loads

Introduction:

Avans DNA Marker is 19 DNA fragment ladder covering wide range of 100-10,000 bp sizes. Conveniently delivered in ready-to-use format and ideal for agarose gel electrophoresis, Northern, Southern, and related DNA analysis applications.

Protocol:

Swirl tube well.

Add ~5µl per load onto agarose gel prior to electrophoresis.

Shipping and Storage Conditions :

Shipped in RT

Store in 4 °C for short term

Store in -20 °C for long term

