

Piggy™ Tube Rotator

The Piggy™ Tube Rotator is designed to work as mixer, tumbler and rotator. Its limitless angles of rotation and variable speed provides user great convenient. Just simply adjusting the holder on main rod, sample can be better mixed. Ideal for molecular hybridization, blood and various samples requiring mixing.

Features

- 10 - 40rpm
- 2 modes selection : variable speed mode & oscillating mode
- 20-99° oscillation angle
- Included 3 kind of rotisseries hold all size of tubes
- Can be used in cold room



PTR-1040-1

PTR-1040-2

PTR-1040-3



PTR-1040

Specifications

Cat. No.	PTR-1040
Speed range	10 ~ 40rpm
Time range	1s - 999min / Continuously
Voltage	AC110 ~ 230V, 50/60Hz
Power	4W
Fuse	250V, 1A, ϕ 5x20
Dimension	W 24 x D 14.5 x H 21 cm
Net weight	1.2kg

Rotisseries

Cat. No.	Description
PTR-1040-1	2 units for 3 x 50ml and 18 x 0.5ml tubes
PTR-1040-2	2 units for 5 x 15ml and 6 x 5ml tubes
PTR-1040-3	2 units for 18 x 1.5/2ml tubes